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User Manual

Pulse Oximeter MATRIX O₂ 20

M. Feingersh & Co.Ltd
2 Dotan St. Yavne, Tel: 08-9999399 The logo consists of the letters "M" and "F" in a stylized, bold font, with a small circular element between them.

INSTRUCTIONS TO USER

Dear users, thank you very much for purchasing the Pulse Oximeter.

This manual describes, in accordance with the pulse oximeter's features and requirements, main structure, functions, specifications, correct methods of transportation, installation, usage, repair, operation, maintenance and storage, etc. As well as the safety procedures to protect both the user and equipment. Please refer to the respective chapters for details.

Please read the User Manual carefully before using this product. The User Manual which describes the operating procedures should be followed strictly. Failure to follow the User Manual may cause measuring abnormality, equipment damage and human injury. The manufacturer is NOT responsible for the safety, reliability and performance issues and any monitoring abnormality, human injury and equipment damage due to user's negligence of operation.

1 SAFETY

1.1 Warnings

A. DO NOT use the oximeter while the patient is being scanned by MRI

B. It is recommended that the oximeter should be inspected before use, when there is obvious damage, stop using the device.

C. Discomfort or pain may appear if using the device ceaselessly. Especially for microcirculation barrier patients, it is recommended that the oximeter should not be used on the same finger more than 30 minutes.

D. For some patients, need a more careful examination of the measurement site, the oximeter shall not be placed in edema or fragile organization

E. The device is not used for neonate and infant.

F. The device is just one of clinical auxiliary equipment, other means of clinical diagnosis and clinical response factors determine the patient's condition must be consolidated under the guidance of a doctor.

G. It is not recommended to use the oximeter in high frequency environment such as electrosurgical equipment.

H. Prevent children from swallowing the oximeter or its accessories children must be accompanied by adult guardian using products.

I. A flexible circuit connects the two parts of the oximeter, DO NOT twist or pull on the connection.

G. Please follow local ordinances and recycling instructions regarding disposal or recycling of the device and batteries.

1.2 CAUTIONS

A. The equipment is designed to measure the percentage of arterial oxygen saturation of functional hemoglobin. Factors that may degrade pulse oximeter performance or affect the accuracy of the measurement include the following:

Do not apply the pulse oximeter on the same arm as a blood pressure cuff, arterial catheter or infusion line(s)(IVs).

- Excessive light, such as sunlight or direct home lighting.
- Moisture in the equipment.
- Finger is outside recommended size range.
- Poor pulse quality.
- Venous pulsations.
- Anemia or low hemoglobin concentrations.
- Cardio green and other intravascular dyes.
- Carboxyhemoglobin.
- Methemoglobin.

- Dysfunctional hemoglobin.
- Artificial nails or nail polish.

B. Please read the measured value when the waveform on screen is equally and steadily-going, the measured value is optimal value, and the waveform at the moment is the standard one.

C. The light (the infrared is invisible) emitted from the equipment is harmful to the eyes, so the user and the maintenance man should not look at the light source.

2 THE BASICS

The pulse oxygen saturation is the percentage of HbO_2 in the total Hb in the blood, so-called the O₂ concentration in the blood. It is an important bio-parameter for respiration. A number of diseases relating to respiratory system may cause the decrease of SpO₂ in the blood, furthermore, some other causes such as the malfunction of human body's self-adjustment, damages during surgery, and the injuries caused by some medical checkup would also lead to the difficulty of oxygen supply in human body. And the corresponding symptoms would appear as a consequence, such as vertigo, impotence, vomit etc. serious symptoms might bring danger to human's life. Therefore, prompt information of patients SpO₂ is of great help for the doctor to discover the potential danger, and is of great importance in the clinical of medical field.

2.1 Principle

Principle of the oximeter is as follows. An experience formula of data process is established taking use of Lambert Beer Law according to Spectrum Absorption Characteristic of Reductive hemoglobin Hb and oxyhemoglobin HbO₂ in glow & near-infrared zones. Operation Principle of the device is Photoelectric Oxyhemoglobin Inspection Technology is adopted in accordance with Capacity Pulse Scanning Recording Technology, So that two beams of different wavelength of lights can be focused onto human nail tip through perspective clamp finger-type sensor. Then measured signal can be obtained by a photosensitive element information acquired through which will be shown on screen through treatment in electronic circuits and microprocessor.


Figure 1 Oximeter schematic diagram

2.2 INTRODUCTION

2.2.1 Intended Use

The Pulse Oximeter is intended to measure functional arterial oxygen saturation (SpO₂) and pulse rate of adult and children patients.

WARNING:

This pulse oximeter is intended for use only by clinical professionals or under their guidance. It must only be used by persons who have received adequate training in its use. Anyone unauthorized or untrained must not perform any operation on it.

CAUTIONS:

- This pulse oximeter is intended for use in hospital clinical institution healthcare community.
- The pulse oximeter is NOT designed for newborn and infant. For adults and children, it is recommended that the finger thickness should be between 8-25.4mm.

2.2.2 Components

The oximeter consists of probe, electronic circuits, and display and plastic enclosures.

NOTES

- The probe is the hole in the middle of the equipment to which the finger insert.
- The probe is the Applied Part of the equipment.

2.2.3 Features

- The pulse oximeter is small in volumes, light in weight and easy to carry.
- Low in power consumption, 600 spot-checks on two AAA batteries.
- One button and easy to operate.
- Auto turn off within 8 seconds when there is no signal.

2.3 Front View


Figure 2. Front View of FS10A


Figure 3. Front View of FS20A

2.4 Functions

Function	FS10A	FS20A
Display	LED	OLED
Spo2 parameter measurement	Y	Y
Pulse rate parameter measurement	Y	Y
Bar graph display	Y	Y
Battery display	Y	Y
Automatic power-off function	Y	Y
Pulse beat sound	Selective	Y
Pulse sound On/Off function	Selective	Y
Pulse waveform display	—	Y
Four direction display	—	Y

2.5 Symbols

Symbol	Definition	Symbol	Definition	Symbol	Definition	Symbol	Definition
%SpO2	The Pulse Oxygen Saturation (%)	[SN]	Serial Number	[]	WEEE(2002/96/EC)	[]	This way up
PR	Pulse Rate (BPM)	[]	The device has no Alarm System	[]	Keep dry	[]	Stacking limit by number
IPX1	The product is protected against harmful effects of dripping water per IEC 60529.	[]	Date of Manufacturer	[]	Fragile, handle with care	[]	Humidity limitation
	Caution, consult accompanying documents	[CE]	This item is compliant with Medical Device Directive 93/42/EEC	[+]	Battery positive electrode	[]	Temperature limitation
[]	BF type applied part	[]	Manufacturer	[—]	Battery cathode electrode	[]	Atmospheric pressure limitation

3 BATTERY INSTALLATION

A. Put the two AAA batteries into battery compartment in correct polarities.

B. Push the battery cover horizontally along the arrow shown as right.

WARNINGS:

Battery polarities should be correctly installed, otherwise, damage may be caused to the equipment.

Please remove the batteries if the equipment will not use for a long time.


Figure4. Battery Installation

4 OPERATING GUIDE

4.1 Application Method

A. Open the battery cover, and put the two AAA batteries into the battery compartment in correct polarities, then replace the cover ;

B. Gently press the bottom of the equipment and open the probe, then insert one finger into the probe;

C. Press the button to turn the equipment on, and the measure interface will appear 3 seconds later ;

D. After about 10 seconds, the measurement result can be read directly from the display screen;

E. Before reading the parameters, make sure that stable waveform of the pulse oximeter interface has sustained more than 4 second;

F. * After turning the equipment on, the pulse sound is open. Long press the button for 2 seconds can close the pulse sound. Long press the button can open the pulse sound again;

G. * In the normal measure state, short press(less than 2 seconds) the button can switch the display direction of the screen;

H. The equipment can be turned off automatically within 8 seconds when the finger left the probe.

Notes:

- The pulse sound are optional function.
- Only FS20A can switch the display directions.


Figure 5. Operation Guide

4.2 Attention for Operation

A. Before use check and confirm that the people or finger size were applicable;

B. Before use check and confirm that the environment should be non-combustible material, as well as to avoid high or low temperature and humidity, but also need to pay attention to the following:

a) To avoid glare and direct sunlight exposure;

b) To avoid radiation infrared or ultraviolet radiation;

c) Avoid contact with the organic solvent, mist, dust, corrosive gases;

C. The finger should be in a proper position, the patient should insert the finger into the probe deeply

- D. The equipment should not be used at a location or limb tied with arterial canal or blood pressure cu or receiving intravenous injection;
- E. The equipment may not work normally on microcirculation barrier patients Warm or rub the finger, or re-position the equipment could improve the measurement.
- F. Don't shake the nger and tr y to keep the patient still during the measurement;
- G. The ray between photo detector and light emitting diode should across patient's arteriole;
- H. The patient should not use enamel or other makeup;
- I. Avoid to insert a wet nger into the probe;
- J. A exible circuit connects the two halves. Do not drag or pull the exible circuit or overextend the equipment's spring;
- K. Do not hang the lanyard from the equipment's exible circuit;
- L. Stop using the equipment if the liquid ows into the inside of the equipment, and wait for the machine to dry.
- M. The pulse waveform is un-normalized waveform, which can be used for signal quality indicator indirectly. The parameter readings may not reasonable when observed disturbed pulse waveform from the pulse oximeter interface.

5 SPECIFICATIONS

5.1 Classication

5.4 Environmental Specications

Type of protection against electric shock:	Internally powered equipment)
Degree of protection against electric shock	Type BF-Applied part (non-debrillation proof)
Operat ng mode	Spot checking
Degree of protection against hazards of explosion	Ordinary equipment: Not protected
Equipment type	Fingertip oximeter
5.2 Measurement Specications	
Spo2 declared accuracy	
Range (σ^*)	7 0%~1 0 %: \pm 2 digits
Resolution	0% ~ 69%: unspecified
PR declared accuracy	1%
Range (σ^*)	2 5~2 5 0 : \pm 3 digits
Resolution	1 bpm

5.3 Power Requirements

Specification of alkaline batteries	Two AAA
Operating current	Less than 3 0mA
Run time	600 spo t checks on two full power batteries at ambient temperature 25

5.4 Environmental Specications

Temperature

Operating	+41°to+104°F/5°to + 4 0°C
Storage/Transportation	- 40°to +140°F/-40° to + 60°C

Humidity

Operating	10~ 95%, noncondensing
Storage/Transportation	10~ 95%, noncondensing

Atmosphere Pressure

Operating	70~106 kpa
Storage/Transportation	50~107.4 kpa

5.5 Physical Specifications

Width*Height*Depth

62 35 31 mm

Weight

60 g (including the batteries)

5.6 Display

	FS10A	FS20A
Display type	LED	Dual_color OLED, 0.96", 128 64 pixel
Display content	SpO2%, Pulse rate, Battery indicator, Bar graph	Pulse waveform, Heart beat mark

Notes:

- σ represents approximately 68% of measurements.
- Specifications are subject to change without noticed.

6 MAINTENANCE,CLEANING,DISINFECTION

6 Maintenance, Cleaning, Disinfection

6.1 Maintenance

The equipment's design life expectancy is about 2 years, keep your equipment and accessories free of dust and dirt, and follow these rules:

- A. Please clean the equipment before use according to chapter 6.2; Remove the batteries inside the battery cassette if the equipment will not be operated for a long time;
- B. Replace the batteries in time when the battery voltage indicate lamps were empty;
- C. It is recommended that the equipment should be kept in a dry environment with no corrosive gases and good ventilation anytime. The moisture and high-light environments will affect its lifetime and even might damage the equipment.
- D. It is best to preserve the product in a place where the temperature is between -20 to 60° and the relative humidity is less than 95%.
- E. The packed equipment can be transported by ordinary conveyance. The equipment not be transported mixed with toxic, harmful, corrosive materials.

WARNING

No modification of this equipment is allowed.

6.2 Cleaning

Your equipment should be cleaned on a regular basis. If there is heavy pollution or lots of dust and sand in your place, the equipment should be cleaned more frequently. Before cleaning the equipment, consult your hospital's regulations for cleaning the equipment. Recommended cleaning agents are:

- a) Mild soap (diluted).
- b) Sodium hypochlorite bleach (diluted).
- c) Hydrogen peroxide (3%).
- d) Ethanol (70%).
- e) Isopropanol (70%).

To clean your equipment, follow these rules:

- a) Shut down the pulse oximeter ;
- b) Clean the display screen using a soft, clean cloth dampened with a glass cleaner ;
- c) Clean the exterior surface of the equipment and probe using a soft cloth dampened with the cleaner ;
- d) Wipe off all the cleaning solution with a dry cloth after cleaning if necessary;
- e) Dry your equipment in a ventilated, cool place.

To avoid damage to the equipment, follow these rules:

CAUTIONS

- Always dilute according the manufacturer's instructions or use lowest possible concentration.
- Do not immerse part of the equipment in the liquid.
- Do not pour liquid onto the equipment or accessories.
- Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).
- If you spill liquid onto the equipment, contact us or your service personnel.

6.3 Disinfection

Disinfection may cause damage to the equipment and is therefore not recommended for this pulse oximeter unless otherwise indicated in your hospital's servicing schedule. Clean the pulse oximeter before disinfecting it.

The recommended disinfectants includes: ethanol 70%, isopropanol 70%, glutaraldehyde-type 2% liquid disinfectants.

CAUTION

Never use ETO or formaldehyde for disinfection.

6.4 Disposal

Dispose of the pulse oximeter in accordance with local environment and waste disposal laws and regulations.

7 ACCESSORIES

One lanyard.

Two AAA batteries.

One Little bag.

One user manual.

One certificate card.

Note:

For particular configuration of accessories please refer to the product package list.

TROUBLESHOOTING

WARNINGS

- Necessary maintenance must be performed by qualified service personal ONLY.
- Users are NOT permitted to maintain the equipment by themselves.
- There are NO replaceable components in the equipment.

Trouble	Possible Reason	Solution
The device can not be turned on.	The battery is drained away or almost drained away.	Please replace battery.
	The battery installation is incorrect.	Install the battery over again.
	The malfunction of the device.	Please contact the local service center.
The display is off suddenly.	The device is set to shut down automatically in 8 seconds when there is no correct physiological signals	Normal
	The battery is almost drained away.	Please replace battery.
Pulse sound can't be turned off	The key is bad.	Check the key and press again.
	The press time is not right.	Make sure the press time is 2~3 seconds.
The Spo2 and Pulse Rate are not displayed stably.	The luminescent or photoelectric tube is sheltered by some object.	Check the luminescent and photoelectric tube.
	The finger is shaking or the patient is moving.	Let the patient keep still.
	The finger is not placed inside deep enough.	Place the finger properly and try again.
	The finger's size is too big or too small.	Select the correct size finger to measure.
	Excessive ambient light	Avoid the excessive ambient light irradiation.
	Pulse rate value of the cyclical fluctuations.	The measurement is normal, and the patient is arrhythmia.
The Spo2 and PR are not displayed normally	The finger is not properly position	Place the finger properly and try again.
	The patient's SpO2 is too low to be detected	Try again, Go to a hospital for a diagnosis if you are sure the device works all right

9 ELECTROMAGNETIC INTERFERENCE

The sharp increase in the number of radio transmitting apparatus or other sources of electrical noise in healthcare environment, these devices are too close or the transmit power is too strong, causing serious interference may interrupt the pulse oximeter to work normally. Which may cause electromagnetic interference sources include:

- Electronic surgical instruments
- Mobile telephone
- Vehicle radio communications equipment
- HDTV

Electromagnetic interference may cause the pulse waveform confusion, large deviation from the measured value, the reading is unstable, or other errors. If this happens, the user should identify sources of interference, and to take the following measures to eliminate interference:

- Turn off nearby devices and then turn on to find the interference device.
- Change the direction or position of the interference device.
- Increase the distance of the interference device and oximeter ;

APPENDIX A

The equipment complies with the requirement of standard EN 60601-1-2:2007 "Electromagnetic Compatibility - Medical Electrical Equipment".

Guidance and manufacturer's declaration – electromagnetic immunity

The FS10A/FS20A Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the FS10A/FS20A Pulse Oximeter should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment –guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity

The FS10A/FS20A Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the FS10A/FS20A Pulse Oximeter should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the FS10A/FS20A Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. $d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P} \quad 80\text{MHz to } 800\text{ MHz}$ $d = \left[\frac{7}{E_1} \right] \sqrt{P} \quad 800\text{MHz to } 2.5\text{ GHz}$ Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic 

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – electromagnetic emissions

The FS10A/FS20A Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the FS10A/FS20A Pulse Oximeter should ensure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The FS10A/FS20A Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	N/A	The FS10A/FS20A Pulse Oximeter is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations flicker emissions IEC 61000-3-3	N/A	

Recommended separation distances between portable and mobile RF communications equipment
and the Medical FS10A/FS20A PULSE OXIMETER

The FS10A/FS20A Pulse Oximeter is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Medical FS10A/FS20A PULSE OXIMETER can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the FS10A/FS20A Pulse Oximeter as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter(W)	Separation distance according to frequency of transmitter(m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
0,01	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$ /	0.02	0.03
0,1	/	0.06	0.11
1	/	0.18	0.35
10	/	0.57	1.1
100	/	1.8	3.5

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.