
Manual do Usuário User Guide

www.wtavet.com.br

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ALLICE Embryo Freezer – WTA

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ADVERTÊNCIA

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Index

WARNING	18
Copyright	18
PRESENTATION.....	19
GENERAL DESCRIPTION	19
ACCESSORIES.....	19
TECHNICAL CHARACTERISTICS.....	20
FUNCTION SPECIFICATION.....	21
IMPORTANT PRECAUTIONS	22
Repairs	22
Cares and Transportation.....	22
Cleaning	22
Falling Objects Inside the Equipment	22
Power Cable.....	22
Rain and Thunderstorms.....	22
SETTING, CONFIGURATION AND OPERATION	23
EMBRYO STRAWS ACCOMMODATION.....	26
NOTICES.....	26
WARRANTY	27
PROBLEM RESOLUTION	30
TECHNICAL ASSISTANCE SERVICE/CUSTOMER SERVICE	30

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USER GUIDE

WARNING

To prevent accidents and/or damage, carefully and completely read and follow this instruction guide before using the product.

Warning for dangerous
voltage presence.

Warning for important
operating instructions.

Copyright

ALLice - Embryo Freezer - user guide

WTA, Cravinhos, State of São Paulo, Brazil.

The content of this guide is protected by copyright laws and cannot be reproduced, distributed, disclosed nor stored without prior written permission of WTA - Watanabe Tecnologia Aplicada Eireli - EPP.

In order to obtain authorization, additional information or clarification, please contact our Technical Assistance Department by e-mail, letter or phone.

Phone: +1 (979) 676-1968

e-mail: sales@wtavet.com.br

*Opening hours: From Monday to Friday, 8:30 am to 5:30 pm (GMT -03:00)

©WTA - Watanabe Tecnologia Aplicada Eireli - EPP - All rights reserved.

Comments

The information, illustrations and specifications contained in this user guide assume the latest features of the product available at the time of its publication.

Due to the use of new technologies and raw materials, WTA reserves the right to change the design and/or technical specifications of the product always as necessary and without prior notice, without incurring in any obligations.

This manual is an integral part of the product. Therefore, the consumer should save it for any queries.

PRESENTATION

Congratulations on your choice. You have just purchased a WTA product, a brand that adds quality and technology in what it manufactures.

WTA is a Brazilian company specialized in solutions for animal reproduction. Our products are developed for working with cattle, sheep and horses, with recognized approval by the most renowned specialists in the world.

Its anodized aluminum frame offers, besides an elegant design, conditions that facilitate its internal and external cleaning.

For the best performance of the product we recommend reading this manual.

GENERAL DESCRIPTION

The ALLice is an embryo freezer that has a micro controlled system that can accurately control the required temperature variations in embryo freezing processes.

It has colorful touchscreen display, intuitive menus and easy operation.

The WTA Embryo Freezer has an internal battery with no-break function; Temperature control with high precision and 3 freezing protocols that can be adjusted.

ACCESSORIES

When purchasing your equipment make sure that your accessories are intact and included in the package:

1 control module

1 embryos block
freezer

straws handle

Power adapter

Obs. Liquid Nitrogen reservoir not supplied.




Note: If any accessory is missing or damaged, contact with the WTA customer service department, within three (3) days from the date of delivery of the product.

TECHNICAL CHARACTERISTICS

- Aluminium manufactured
- Control unit dimensions: 260 x 200 x 105mm.
- Freezing block dimensions: (D x H) 85 x 150mm.
- Capacity: 18 straws of 0.5 mL or 36 straws of 0.25 mL (2 straws / well).
- Approximate gross weight: 6,84lb.
- Maximum power: 64,5 W.
- Power AC supply: 100~240Vac - 50/60 Hz.
- Equipment DC supply: 15 Vdc – 4,3A.
- Working temperature: -58,0 to 86,0°F (-50,0 to +30.0°C).
- Temperature variation: $\pm 0,36^{\circ}\text{F}$ ($0,2^{\circ}\text{C}$).
- Battery life: 3 hours after full charge.
- Product limited warranty: 2 years.

FUNCTION SPECIFICATION

1. Switch ON/OFF
2. Display
3. Straws handle
4. Embryo freezing block (rack)
5. DC power inlet
6. USB inlet
7. Block resistance supply cable (4 pins)
8. Temperature sensor cable (5 pins)
9. Led indicator
10. Nitrogen liquid level

IMPORTANT PRECAUTIONS

Repairs:

Do not open the device. All repairs and maintenance should be performed by authorized professionals or by WTA service.

Cares and Transport:

Avoid shaking during transport which can damage it. Make sure that the equipment is properly seated and protected against any accidents.

Cleaning:

Clean the exterior with a soft cloth slightly moistened with water. Do not use alcohol, thinner or other solvents as they will cause damage device's finishing.

Falling Objects Inside the Equipment:

Do not place recipients or cups containing liquid, or coins, paper clips or other metallic objects on the control module. If any liquid or any objects falls inside the equipment, turn off the power immediately and contact the service center.

Power Cable:

Do not step on the power cable or put any kind of objects (furnishing, etc.) over it.

Rain and Thunderstorms:

In storms with electric discharge events disconnect the equipment from the AC power.

SETTING, CONFIGURATION AND OPERATION

Install the equipment over a leveled and dry even surface.

Connecting the Equipment

For safety, first connect the output DC power into the control module input (05) and then connect the power supply to the 100 ~ 240Vac mains.

Turn ON / OFF

To turn the equipment on, press the switch (01) to the ON position and to turn it off press again to the OFF position.

The "LED" Indicator (9) allows to monitor the charging of the battery:


- flashing red = external source (equipment without battery);
- solid red = charging battery;
- solid green = fully charged battery.

The display icon  (2) allows you to monitor the battery charging progress.

When the equipment is running on battery power the "LED" will be off.

Configuring the Equipment

After the animation presentation the "Home" screen is displayed.

A simple icons touch triggers the following screen, in which the saved information can be viewed or made new adjustments.

Displays information: - Owner

- Firmware versio

- Serial number

- Date of Manufacture

Exclusive WTA Service use, together with USB input (6).

Setting date and time:

Touch the window and enter date and time values in the format indicated. Press  to confirm.


Settings and Execute Freeze Curve


When touch the "Run Curve" icon a new screen will be displayed with three freezing protocols, identified by buttons **1**, **2** and **3**, which can be adjusted according to the user's needs.

After selecting the desired protocol, a new screen will be displayed with all information

To change the parameters of the freeze curve, touch inside the windows. The adjustment screens will open. Touch "+" or "-" to adjust the value. Touch  to confirm.

Adjust the other parameters:

- (1) - Target 1; Rate 1; Hold Before Seeding; Hold After Seeding;
- (2) - Target 2; Rate 2; Hold Time 2;
- (3) - Target 3; Rate 3; Hold Time 3.


Note: The ALLice is designed with 3 freezing protocols. Changes of values will be saved and used again when the equipment is turned on. **Recommended to write down the values before changing them.**

After the freezing setting is completed, fill the tank with liquid N₂ and touch **START**.

Place the freezing block (4) inside the tank with liquid Nitrogen and avoid exceeding the set level (10). Freeze block cables (7 and 8) must not be interlaced. During freezing the cables will be frozen due to Nitrogen vapor. Avoid touching or moving them at risk of rupturing.

Touch OK to start freezing.

EMBRYOS STRAWS ACCOMODATION

Follow the freezing block temperature through the display (2). When the start temperature is reached, the insert straws message will be displayed.

Insert the straws, and touch OK.

The programming will go to target 1 with rate 1.

On reaching target 1, the pre-seeding time message will be displayed.

After completing the pre-seeding time the message to start seeding will be displayed.

The programming will be waiting to finish the manual seeding.

After confirming the manual seeding, the post seeding will be started.

After completing the "seeding" manual, start the post seeding.

After the post "seeding" the programming will go to target 2 with rate 2.

Arriving on target 2 the temperature will be maintained during stabilization 2.

After stabilization 2 is complete, the programming will go to target 3 with rate 3.


Arriving on target 3 the temperature will be maintained during the stabilization time 3.

At the end of stabilization 3 the finalization message will be displayed.



Touching OK will disable the temperature control.

NOTICES

During the freezing procedure the display will show status messages, for example:

Warning message will be displayed if any cable is disconnected.

WARRANTY

WTA - Watanabe Applied Technology Eireli - EPP follows strict internal quality control standards in the products manufacture. In addition to the legal guarantee of 90 (ninety) days, provided for in the Consumer Protection Code, it guarantees its customers an additional 21-month Warranty, counted from the date of issue of the purchase invoice, provided that the equipment has been installed and according to the guidelines contained in the User's Manual.

Even if ownership of the equipment is transferred, the warranty will continue to be valid according to the terms and conditions specified above.

This product has been carefully evaluated at every stage of the manufacturing process. This warranty will always accompany you, valid only in the country where the product was purchased.

It is understood by Additional Warranty, the free repair of the product and / or the replacement of parts that, according to technical opinion of the WTA, present a manufacturing defect.

Only WTA Technical Assistance professionals are entitled to repair defects for the guarantee, by presenting the sales tax (original) to the first buyer / consumer.

The Warranty Service

- Assistance in WTA technical assistance.
- Parts used for maintenance, subject to the exceptions provided in the item "Exclusion of Warranty".

The product should be sent directly to the WTA observing that:

- The product must be properly packed, if possible in its original packaging, together with the manual, cables and other accessories that have been delivered at the time of purchase.
- The product will be pre-evaluated by WTA technicians prior to repair.
- Shipping, freight and packaging costs will only be charged to WTA during the legal guarantee period (90 days). After this period such costs will be the customer's responsibility.

Voiding the Warranty

This warranty does not cover defects or defects arising from the following situations:

- Installation and use of the product in an improper way or in disagreement with the respective guidelines and information described in this Manual, including installation, external adjustment, cleaning and use of accessories not in accordance with technical specifications.
- Defects caused by accidents, falls, atmospheric phenomena (rains, floods, electric discharges), excessive humidity (water, liquids in general), salinity, sinister (theft or theft), loss or improper transport of the product.
- Installation, adjustment, repair or servicing by persons outside WTA service.
- Invoice and / or serial number adulterated or withdrawn.
- Modifications, alterations or violation of the original product.
- Product connected to the mains improper or subject to excessive fluctuations.
- The product is purchased in balance, showcase or demonstration, this being characterized by means of identification in the invoice of purchase and / or price lower than the one practiced in the market.

Important: WTA is not liable for damage to equipment during transportation. When purchasing your product make sure on the delivery act that it is not defective or with apparent defects. If you notice any irregularities, return the product immediately.

This Warranty Statement Does Not Include

- Periodic reviews.
- Displacement and expense of technician by authorized WTA service.

Responsibility

The WTA has no control or influence over the conditions under which the equipment is used, as well as the methodologies employed. Thus, WTA assumes no responsibility for the results, use and / or performance of the equipment (including: damages for loss of profits, interruption of business activity, loss of business information or other pecuniary losses). WTA expects that the equipment will be used by trained personnel. The purchaser acknowledges that the equipment is subject to all terms and conditions described in this manual.

PROBLEM RESOLUTION

Occasionally many apparent defects may be something simple to solve. If you notice any problems with the equipment, check the following procedures.

Equipment Is Turning Off When Starting a Rate

- Check the battery level.
- Check that the power cable is connected correctly.
- Turn the machine off, wait 5 seconds, then on again.

Freeze Curve Is Not Operational

- Check that the cables are properly connected.
- Check the level of liquid Nitrogen in the reservoir.

If the problem persists, contact your service representative.

TECHNICAL ASSISTANCE SERVICE/CUSTOMER SERVICE

If there is any doubt or if it is necessary additional information about warranty, assembly, installation, adjustment or utilization of the product, please contact our Customer

Service Department:

Opening hours: From Monday to Friday, 8:30 am to 5:30 pm (GMT -03:00)

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Todos os direitos reservados | All rights reserved