



Analytical and Biochemical Method Development Examples

Bio-Technical Resources (BTR)

www.biotechresources.com

1035 S 7th Street, Manitowoc, WI 54220

Phone: (920) 684-5518


Fax: (920) 684-5519

Email: info@biotechresources.com

Example 1: Optimization of Process Conditions for Linoleic Acid Bioconversion

- ❖ Linoleic acid C9 isomerase converts linoleic acid (*c9,c12*, 18:2) to conjugated linoleic acid single isomer (*t10,c12*, 18:2), a high value dietary supplement
- ❖ A recombinant LAC9 isomerase was produced in *E. coli* by fed-batch fermentation to prepare enzyme catalyst for the development of a bioconversion process
- ❖ Bioconversion optimization:
 - Increase conversion efficiency
 - Improve fatty acid mass balance

Fig. 1-1 Unoptimized Bioconversion Process




LAC9 isomerase sample was incubated with linoleic acid (LA) and FAMES were extracted and analyzed by GC. Conjugated linoleic acid (CLA) was formed rapidly within 3-4 hr.

- CLA was unstable during the process: 70% CLA produced earlier was lost by 16 hr.
- Poor mass balance: total level (LA+CLA) dropped to 15% within 16 hr.



Fig. 1-2 Improved Bioconversion Process



Recombinant *E. coli* cells expressing the LAC9 isomerase or partially purified enzyme were incubated overnight with linoleic acid (LA) in reaction Buffers A and B to produce conjugated linoleic acid (CLA). FAMES were extracted and analyzed by GC.

- ❖ Optimized conditions (Buffer B) lead to much higher levels of CLA production and mass balance than Buffer A in both whole cell biotransformation and enzymatic conversion.
- ❖ Mass balance of bioconversion reached ~100% in Buffer B using purified enzyme.

BTR



Example 2. Assay Optimization to Reduce Variability of Bioluminescence Detection

❖ Client needs


- Develop a proprietary process to produce a bioluminescent protein by recombinant *E. coli* fermentation for pharmaceutical applications
- Project scope: vector, strain and fermentation process development

❖ Analytical support

- Protein quantification by bioluminescence, SDS-PAGE and HPLC
- Challenge: bioluminescence assay is very sensitive but suffers from high variability (up to 30%)
- Assay optimization: variability reduced to within 5%



Fig. 2-1 Optimized Bioluminescence Assay



Dilution	Enzyme (ng/ul)	RLU/s	stdev (n=4)	RSD (%)
1/2	2.0	1,976,019	27,482	1.4%
1/4	1.0	1,066,837	48,520	4.5%
1/8	0.5	576,985	23,955	4.2%
1/16	0.3	308,003	11,660	3.8%
1/32	0.1	154,239	4,226	2.7%
1/64	0.1	76,841	1,099	1.4%



Example 3: Development of a HPLC Method to Monitor Protein Production in Fermentation

❖ Client needs

- Fermentation process development for expression and secretion of a protein using *Bacillus subtilis*
- Analytical need: a reliable method for product quantification

❖ Method Development

- Client unoptimized method had poor reproducibility: peak split and peak area reduction after ~10 injections using guard column
- Problems overcome by cleaning samples using Bio-Spin columns or 96-well microdialysis plates: laborious and added cost
- A new HPLC protocol was developed by changing elution profile

❖ Outcome

- No need for laborious sample cleaning
- Results are highly reproducible



Fig. 3-1: Protein Standard Peak Tailing and Splitting After Multiple Sample Injection (Old Method)






Fig. 3-2: Samples Cleaning Using Bio-Spin Column



Chromatograms show the peak of Protein Standard (0.5 g/L) injected before injection of any fermentation samples (A), peak split after 12 injections of untreated samples (B) and no peak split after 119 injections of samples cleaned using Bio-Spin column (C)




Fig. 3-3 Fermentation Samples Spiked with Protein Standard



A fermentation sample (A) and the sample spiked with protein standard (B) show sharp peaks on the chromatograms using the improved HPLC method



Fig. 3-4 Performance of the New HPLC Method



- ❖ Linear range of detection 0-1.0 g/L
- ❖ Consistent results after large numbers of sample injection