


## Introduction

Influenza viruses are responsible for seasonal epidemics each year in the United States that result in an average of **200,000** hospitalizations and tens of thousands of deaths.

**Accurate and early diagnosis of influenza viral infections is critical.**

Using simple and efficient pre-concentration workflows combined with existing molecular test formats can improve detection at earlier time points.



<https://molecular.diasorin.com/>

## Nanotrap® Particle Technology

**Nanotrap® particles are engineered hydrogel particles that use chemical affinity baits to capture, concentrate, and preserve low abundance analytes from biological samples.**

### Core:

- Highly porous hydrogel facilitates rapid analyte exchange
- Functionalized with affinity bait to capture analytes of interest

### Affinity Bait:



- Selected to have very high affinity for a specific class of analyte, such as proteins, peptides, nucleic acids, exosomes, viral particles, or bacteria

### Shell (optional):

- Second layer of hydrogel performs size exclusion.
- Can be chemically modified for additional functionalization.

### Molecular Functionalization of Shell (optional):

- Selected to provide additional functionality, such as to repel charged analytes



## Nanotrap® particles for viral enrichment

Nanotrap® particles concentrate viral pathogens in biological fluids for improved detection by downstream diagnostic assays such as enzyme-linked immunosorbent assay (ELISA), next generation sequencing (NGS), lateral flow assays (LFA), plaque infectivity assays, and quantitative reverse-transcription polymerase chain reaction (qRT-PCR).

**Table 1. Nanotrap® Particles for Flu A Viral Capture**

NT Catalog No.	Product Description
CN1030	Reactive Red 120 Core Hydrogel Particles
CN3080	Magnetic Reactive Red 120 Core Hydrogel Particles

P. Andersen<sup>1</sup>, D. Goldfarb<sup>1</sup>, J. Warsh<sup>1</sup>,  
D. Munns<sup>1</sup>, E. Porter<sup>1</sup>, B. Lepene<sup>1</sup>, R. Barbero<sup>1</sup>  
<sup>1</sup>Ceres Nanosciences, Inc., Manassas, VA, USA

## Experimental Overview

### Objectives


- To demonstrate that magnetic Nanotrap® particles can significantly improve sensitivity of a flu A molecular assay by enriching flu virus from a larger sample.
- To demonstrate the use of Nanotrap® particles to pre-concentrate Flu A from large volume viral transport media samples instead of a viral RNA extraction column for flu molecular assays.

## Methods Comparison

### QIAGEN Column-Based Method

150 µL of viral transport media spiked with Flu A positive control is mixed with RLT buffer and processed with QIAamp Viral RNA Mini Kit (QIAGEN) and then analyzed with a Genesig Flu A assay.

### Magnetic Nanotrap® Particle-Based Methods



Nanotrap® particle method outperforms QIAGEN columns for extraction of Flu A RNA from spiked samples




Figure 1. QIAGEN viral RNA mini Kit and Ceres Method 1 for H1N1 virus RT-qPCR applications

## Large Volume Sample Enrichment

**Ceres Method 1** enriches virus from large volumes of samples, which brings the low concentration sample into the detectable range, and results in a **19-fold increase** in the detected levels of Flu A at the high concentration.




Figure 2A. Enrichment comparison between QIAGEN columns and using Ceres Method 1 with a larger starting sample volume (2 mL).




Figure 2B. Method comparison between QIAGEN columns alone and Ceres Method 2 with increased sample volume (2 mL).

Combining Nanotrap® particles with QIAGEN columns in **Ceres Method 2** further improves Flu A enrichment resulting in **≥ 8-fold** and **≥ 3-fold** increase in detected levels of Flu A.

## Nanotrap® Particles & Automation

Magnetic Nanotrap® particle-enabled large-volume sample enrichment in **Ceres Method 3** offers streamlined workflows that can facilitate high-throughput automated sample processing on liquid handling platforms.




Figure 3. Ceres Method 3 performance across high and low concentrations of spiked transport medium.

## Conclusions & Future Directions

- Nanotrap® particles can improve flu virus detection by pre-concentrating flu virus from viral transport media samples for RT-PCR detection.
- Methods are amenable to automation.
- Nanotrap® particles preserve pathogens at elevated temperatures. Incorporating Nanotrap® particles in the transport media could further improve sample processing workflows and downstream assay.