

mooyee
Relaxer

Instructions

SAFETY NOTES

Contraindicated as follows

1. People with pacemaker, defibrillator, or other implanted electronic devices.
2. Severe cardiac arrhythmia, angina or heart failure.
3. In pregnancy, children under 12 years old.
4. Pets.

Use under medical supervision

1. On treatment or medication, with severe diseases.
2. With fever, in menstrual cycle or in breastfeeding.
3. Acute/ purulent inflammatory, bleeding tendency.
4. Blood pressure disorder, infectious or skin diseases.
5. Skin allergy.

Attentions

1. Keep out of the reach of children.
 2. The old and those who cannot take care of themselves should be accompanied.
 3. Do not put the Gel Pad in your mouth.
 4. No strenuous exercise while in Physiotherapy.
 5. Keep your skin dry and no wounds on the covered area.
 6. Keep away from metal objects such as necklace or bracelet when using.
 7. Stop using it if you feel uncomfortable.
 8. Power adapter is NOT included in the package.
- Smartphone adapters and computer USB ports are suitable; inappropriate adapters may damage the device.


Package Includes

- Mooyee Smart Massager *1
- Gel Pad *2 pairs
- Storage Bag *1
- USB cable for Charging *1
- Instructions *1
- Warranty card *1

Specification

Material: Non-toxic flexible Silica Gel
 Size: 5.5" * 2.2" * 0.5" (141mm*55mm*12mm)
 Weight: 1.73oz (49.2g)
 Connection: Bluetooth (4.0 or above)
 Battery Voltage: 3.7V
 Battery Capacity: 250mAh
 Charge Voltage: DC 5V/1A
 Operation Temperature: 32°F~95°F (0°C~35°C)
 Duration: 5 hours
 Full Charge Time: 2 hours

Introduction to product keys



① : Power button/ Indicator ② : USB charge port ③ : Reset button
 ④&⑤ : Gel Pad contact point


Normal Use Indicator Status:
 Blue: Standby
 Light off: Power off
 Breathing: In operation

Charging Indicator Status:
 Red: Charging
 Blue: Fully charged

Directions for use


Step One

You can find the application with the name 【mooyee】 in GooglePlay or AppStore.



Step Two

Insert the Gel Pad with the button lock;
 Tap and hold the circle button in the middle for 3 seconds till the device is on;
 Peel off the reusable protector film on Gel Pad;
 Put the device where you need to massage.




Step Three


Open App and turn Bluetooth on;
 Auto Device connection check with Bluetooth;
 Choose your Mooyee device and confirm.

Step Four


Choose your preferred scenario and begin relaxing.



Best Areas



Prohibited Areas



Maintenance

1. Turn off the power after use.
2. Get it fully charged before long-term storage and charge it every two or three months to extend the life of built-in battery.
3. Put the device in storage bag if you don't use it.
4. Use baby wipes to clean the stain and dust on device.
5. Place the protector films on Gel Pads after use to keep them away from stain and dust.
6. Use a small amount of purified water drops to clean the Gel Pad in order to maintain the viscosity of the Gel Pad.
7. You need to replace the Gel Pads with new pair when the viscosity goes away after approximately 50 cycle times of use.
8. Please DO NOT USE or STORAGE this device under the following circumstances:
 - In the sun
 - Dusty
 - In case of fire
 - Hot and Humid
 - Electromagnetic Environment

Basic Troubleshooting Guide

Fail to turn on or off	Charge the device. If nothing happens, use a needle-like object to press the reset button in the hole. If the problem still can't be solved, please contact our customer service team.
Device not found	Make sure the Bluetooth and the device is on. If they still can't be connected, please do as follows: 1. Close the Smartphone app; 2. Turn off Bluetooth; 3. Restart the device; 4. Turn Bluetooth on, open application; 5. Search the device. 6. If it can't be found, please try steps in order 1-5 a few more times. If the steps above still can't solve the problem, please contact our customer service team.
Automatically disconnected	Possible causes: 1. Too long distance or barriers between the Smartphone and the device. 2. Insufficient battery power. 3. High-power Bluetooth devices working around. You need to: 1. Keep it within 10 meters, no walls, no large metal barriers between Smartphone and the device. 2. Make sure the battery power is above 20%. 3. Turn off or keep away from high-power Bluetooth devices working around.
Charge failure	1. Please make sure the power adapter output voltage is DC 5V/1A and is working well. 2. Replace the charge cable (or use Smartphone Micro USB cable). If it still can't be charged, please contact our customer services. Notes: In charging mode, the indicator is red or blue (fully charged).
Sting feeling during massage	Attach the Gel Pads to your flat skin part. If the viscosity of Gel Pad fades away, use purified water to clean it. The sting feeling on each proper using occasion may be caused by allergic reaction.

Declaration of hazardous substances

Name of the parts	Toxic or hazardous Substances and Elements					
	Pb	Hg	Cd	Cr(VI)	PBB	PBDE
Shell	○	○	○	○	○	○
PCB	○	○	○	○	○	○
Battery	○	○	○	○	○	○
Gel Pad	○	○	○	○	○	○
USB Cable	○	○	○	○	○	○
Package	○	○	○	○	○	○

○: Indicates that this toxic or hazardous substance contained in all of the homogeneous materials for this component is below the limit requirement in GB/T 26572-2011.
 ×: Indicates that this toxic or hazardous substance contained in at least one of the homogeneous materials for this component is above the limit requirement in GB/T 26572-2011.

All components in the table above meet the 65/2011-/EU instruction of EU RoHS Recast.