

Therapy Support Device - Brief Article

Nursing Homes Magazine, [October 2000](#)

The Second Step Therapeutic Ambulator/Gait Harness System completely encircles the resident with a comfortable padded frame. A supportive gait harness attaches within the frame and provides added padding to protect soft tissue. It has handholds on both sides and the back to help therapists manage residents more safely and easily. The system allows for greater control of a resident's mobility, which helps to facilitate additional therapy opportunities. The lower support works in conjunction with a thoracic belt to distribute weight more evenly.

The Therapeutic Gait Harness uses a quick-release buckle, three handholds and a broad interface between the resident and the harness. The lower harness, which connects with a detachable zipper, secures the lower portion of the body with straps around the thighs. The harness is available in small, medium and large sizes.