

A SINGLE AXIS KNEE with HYDRAULIC regulation of stance phase and swing phase, equipped with a manual lock.

- · Strong high-end knee for dynamic transfemoral amputees.
- · Navigate everyday life activities on all types of terrain; including slopes and stairs.
- · Offers smooth ambulation and the ability to safely descend stairs step over step.

Increased patient dynamics

Optional manual lock

Powerful hydraulic brake

Extension assist

Wide braking range

Maximum patient's weight: 330lbs. (150kg)



STANCE: Adjustment of sensitivity for switching to braking mode.

RESISTANCE: Adjustment of braking power.

LOCK: Manual hydraulic lock switch.

FLEXION: Adjustment of flexion angle.

EXTENSION: Dynamic control of swing phase and terminal impact.



Adjustment wrench (included)

TECHNICAL FEATURES

Item Number for Order	1P130	1P130-KD
Weight	2.79lbs. (1265g)	2.83lbs. (1285g)
TH/BH/H (in mm)	232/177/36	222/174/33
Materials	Aluminium, Stainless steel, Titanium	
Knee flexion	120°	
Proximal Connector	Pyramid	M36 screw
Distal Connector	34mm Tube	
Maximum patient's weight	330lbs. (150kg)	
Activity level	K3-K4	
Water resistance	Weather resistant	

RELATED PRODUCTS

Lamination adapters	1K179	1K179
Adapters	1K183•1K207-HI	D 1K207-HD
	1D41-HD•1K190	1D41-HD•1K190
Tubes	1G01-HD 1G01-HD	
RUSH Foot® Coll	ection	1A600: DynaTrek foot



1236 West Southern Avenue #101, Tempe, AZ 85282 • 855.450.7300







Human First