



Rehabilitation & Mobility Equipment



XXL-REHAB
Aids For Bariatric Users



WELCOME TO THE WORLD OF XXL-REHAB



XXL-Rehab has specialised in supplying aids for bariatric users with the purpose of providing independency, comfort and relief.

About 20 years ago we started to import products for this patient group. We soon realized that most products were just reinforced but not redeveloped to meet the specific demands of bariatric patients.

Now we design and develop our own products. Our products are the result of many years of research in the bariatric field, knowledge obtained from congresses around the world as well as hands-on experience from bariatric users.

We have studied different body types and various mobility patterns as well as geometric studies and the law of gravity. Based on a mixture of this knowledge we design products which maintain the possibility of independency for the user.

Our products are strongly engineered and functionally designed and have a weight capacity of 325 kg. The key words for our product development are innovation and considering the needs of the users.

Our product portfolio consists of wheelchairs, shower commodes, shower benches, patient/dinner chairs, wheeled walkers, walking frames, canes, lift, slings as well as both low air loss and alternating mattresses.

Based in Denmark we supply both to the private and public sector all over the world. Our goal is to provide professional advice and high quality products within a short time frame. Our staff members are highly experienced within occupational therapy and have hands-on experience with bariatric patients.

All people are different but we have in common that we wish for quality of life. At XXL-Rehab we treat all people as were they part of our own family. For us the respect for the individual is very important.



Sales • Hire • Service

SHOWER COMMODE

We have aimed at designing the best shower commode on the market. It is designed for easy access, independence, safety and on the caregiver terms.

The inside width of the chair is 61 or 71 cm, but using a 5 cm high PU foam seat and 7 cm free space to the armrests allows the body mass to fall down, instead of fading out to the sides. The impossible "wide inside but small outside" is obtained. That makes the XXL-Rehab Shower commode passing nearly though all modern doorways.

The footrest is always the weak point when it comes to bariatric products. Normally the producer makes a short cut here, by using a standard solution. We did not only design the footrest to carry 100 kg each, we did more!

The footrest is also side/width adjustable, so access to intimate hygiene is possible and people with the fat-tissue placed on the inside of their leg can sit with legs apart (up to 110 cm) and still have foot support. The brand new footrest is even nicely covered with PU foam.

One of the unique new features is a split toilet seat, a long open hole: 44 cm and 22 cm at the widest point. This to enable a correct use of the chair, it eases access for caregivers and avoids incidents where material ends up in places it wasn't meant to end up!

The ability for bariatric users to reach far back to reach the armrest is usually not present. Therefore we have made the armrests with a length of 64 cm and let it stick out 7 cm from the seat base, ready to welcome and support the hands of mobile users. It enables the user to get into the chair by him/herself. Another option is the butterfly armrests. Butterfly armrests add more than 10 cm to the seating width.

The space between backrest and seat has a clearance of 32 cm. This gives space for a possible bulbous gluteal region, but also enables the caregivers to reach for the sling and position the user better.

- Long armrests provide easy access
- Extra seat width with butterfly armrests



Product video
Scan the code and see
the product in use

STRONGLY ENGINEERED AND FUNCTIONALLY DESIGNED EXCLUSIVELY FOR BARIATRIC USERS!

SPECIFICATIONS

Model 1

Model 2

Maximum weight capacity is	325 kg	
Seat width	61 cm	71 cm
Total width	67 cm	77 cm
Seat depth	56 cm	
Seat height	55 cm	
Backrest	32-52 cm	
Seat to footrest	40-42-44-46-48 cm	
Width of footrest	58-110 cm	
Pushbar from floor	100 cm	
Width between armrests	54 cm	
Height of backrest	48 cm	
Height of footrests	40-48 cm	
Distance from floor to seat underside	44.5 cm	
Total length	117 cm	
Total weight	27 kg	30 kg
Wheels (2 wheels with individual locks)	125 mm	
Seat	Pressure moulded PUR foam	
Frame	Powder-lacquered chromium molybdenum steel	
ISO - classification	ISO 09 12 03	

ORDERING

Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Shower commode incl. standard armrests, footrests and commode bucket	61	56	0150-061-000
XXL-Rehab Shower commode incl. standard armrests, footrests and commode bucket	71	56	0150-071-000
XXL-Rehab Shower commode incl. butterfly armrests, footrests and commode bucket	61 + 14	56	0151-061-000
XXL-Rehab Shower commode incl. butterfly armrests, footrests and commode bucket	71 + 14	56	0151-071-000
Accessories			
Butterfly armrests (add 14 cm to seat width)			0150-099-012
Commode bucket			0140-099-001
Footrest right			0150-099-008
Footrest left			0150-099-009
Calf support / amputation support right			0156-099-030
Calf support / amputation support left			0156-099-031
Shower seat			0141-099-010



Calf support / amputation support



Standard armrests



Side/width adjustable footrests



*Butterfly armrests
add 14 cm to seat width*



Soft shower seat



Sales • Hire • Service

SHOWER COMMODE WHEELER

We have aimed at designing the best shower commode wheeler on the market. It is designed for easy access, independence, safety and on the caregiver terms. The 24" rear wheels provide more directional stability and smooth pushing.

The inside width of the chair is 61 or 71 cm, but using a 5 cm high PU foam seat and 7 cm free space to the armrests allows the body mass to fall down, instead of fading out to the sides. The impossible "wide inside but small outside" is obtained.

The footrest is always the weak point when it comes to bariatric products. Normally the producer makes a short cut here, by using a standard solution. We did not only design the footrest to carry 100 kg each, we did more!

The footrest is also side/width adjustable, so access to intimate hygiene is possible and people with the fat-tissue placed on the inside of their leg can sit with legs apart (up to 110 cm) and still have foot support. The brand new footrest is even nicely covered with PU foam.

One of the unique new features is a split toilet seat, a long opening: 44 cm and 22 cm at the widest point. This to enable a correct use of the chair, it eases access for caregivers and avoids incidents where material ends up in places it wasn't meant to end up!

The ability for bariatric users to reach far back to reach the armrest is usually not present. Therefore we have made the armrests with a length of 64 cm and let it stick out 7 cm from the seat base, ready to welcome and support the hands of mobile users. It enables the user to get into the chair by him/herself. Another option is the butterfly armrests. Butterfly armrests add more than 10 cm to the seating width.

The space between backrest and seat has a clearance of 32 cm. This gives space for a possible bulbous gluteal region, but also enables the caregivers to reach for the sling and position the user better.

- **XXL-Rehab Shower commode wheeler with 24" rear wheels is easy to drive for the attendant**



STRONGLY ENGINEERED AND FUNCTIONALLY DESIGNED EXCLUSIVELY FOR BARIATRIC USERS!

SPECIFICATIONS

Model 1

Model 2

Maximum weight capacity is	325 kg	
Seat width	61 cm	71 cm
Total width	87 cm	97 cm
Seat depth	56 cm	
Seat height	55 cm	
Backrest	32-52 cm	
Seat to footrest	40-42-44-46-48 cm	
Width of footrest	58-110 cm	
Pushbar from floor	100 cm	
Width between armrests	54 cm	
Height of backrest	48 cm	
Height of footrests	40-48 cm	
Distance from floor to seat underside	44.5 cm	
Total length	112 cm	
Total weight	27 kg	30 kg
Front wheels	125 mm	
Rear wheels	24"	
Seat	Pressure moulded PUR foam	
Frame	Powder-lacquered chromium molybdenum steel	
ISO - classification	ISO 09 12 03	

ORDERING

Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Shower Commode Wheeler incl. standard armrests, footrests and commode bucket	61	56	0152-061-000
XXL-Rehab Shower Commode Wheeler incl. standard armrests, footrests and commode bucket	71	56	0152-071-000
XXL-Rehab Shower Commode Wheeler incl. butterfly armrests, footrests and commode bucket	61 + 14	56	0153-061-000
XXL-Rehab Shower Commode Wheeler incl. butterfly armrests, footrests and commode bucket	71 + 14	56	0153-071-000
Accessories			
Butterfly armrests (add 14 cm to seat width)			0150-099-012
Commode bucket			0140-099-001
Footrest right			0150-099-008
Footrest left			0150-099-009
Shower seat			0141-099-010
Calf support / amputation support right			0156-099-030
Calf support / amputation support left			0156-099-031



Calf support / amputation support



Standard armrests



Side/width adjustable footrests



*Butterfly armrests
add 14 cm to seat width*



Soft shower seat



Sales • Hire • Service

SHOWER COMMUNE TILT

The XXL-Rehab Shower Commune Electric Tilt has a 32 degree tilt position and a weight capacity of 325 kg Safe Working Load. This shower commune is designed to provide a comfortable alternative showering system in instances where the prone position might cause discomfort or distress to the user. It may also be used in circumstances of limited space providing a bathing solution where other units may not have room for operation.

Making tasks easier and safer

The XXL-Rehab Shower Commune Tilt provides better access for the caregiver in terms of hygiene, skin and tissue care as well as nursing/dressing feet in standing position, allowing the caregiver to work in a safe and ergonomically sound way.

Improving care, efficiency and communication

This shower commune raises the efficiency of care because it enables a single caregiver to perform the whole showering procedure, as well as a wide range of everyday hygiene tasks. The high level of support provided by the shower commune results in streamlined routines that offer considerable time savings compared with traditional methods. It also enables the caregiver and the user to maintain eye contact for better communication.

Suitable for

The XXL-Rehab Shower Commune Tilt is suitable for patients with reduced muscle tone, pressure ulcers, spasm, involuntary movement, amputations or stroke patients.

- Electric tilt up to 32°
- The high backrest and sturdy armrests provide optimum user support and security
- The reclined position gives full access for the caregiver to most parts of the body and is very comfortable for the user



Product video
Scan the code and see
the product in use



With standard footrests

STRONGLY ENGINEERED AND FUNCTIONALLY DESIGNED EXCLUSIVELY FOR BARIATRIC USERS!

SPECIFICATIONS

Maximum weight capacity is	325 kg
Seat depth adjustable	42 / 48 / 54 cm
Total length (lowest position)	140 cm
Total width	74 cm
Seat width	66 cm
Seat height	51-69 cm
Backrest	24-53 cm
Width of footrest	58-110 cm
Height of footrests	40-48 cm
Electric tilt	0-32°
Width between standard armrests	61 cm
Width between butterfly armrests	76 cm
Width toilet opening	Diameter 24 cm, 35 cm long
Distance from floor to seat underside	44.5 cm
Pushbar from floor	100 cm
Seat to footrest	39-42-45-48 cm
Height seat to armrest	22 cm
Weight	39 kg
Folding armrests	No
Swing-away armrests	Yes
Swing-away footrests	Yes
Removable footrests	Yes
Height adjustable footrests	Yes
Headrest	Yes
4 wheels with individual locks	Yes
Control	Via the control panel
Capacity up/down fully charged battery	Approx. 20 times
Turning radius	158 cm
Frame	Powder lacquered chromium molybdenum steel
Wheels diameter	150 mm
Seat / back / armrests	Polyurethane
IP - classification	IP44
ISO - classification	ISO 09 12 03 - 23

ORDERING

Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Shower Commode Tilt incl. standard footrests and charger	66	42 / 48 / 54	0155-066-000
Accessories			
Butterfly armrests (add 14 cm to seat width)			0150-099-012
Headrest adjustable			0156-099-020
Commode seat, 2 part			0150-099-020
Commode bucket			0140-099-001
Footrest right			0150-099-008
Footrest left			0150-099-009
Calf support / amputation support right			0156-099-030
Calf support / amputation support left			0156-099-031



The handset gives the caregiver control over battery-powered functions



Commode bucket and seat



Headrest



Sales • Hire • Service

SHOWER COMMODE RISE N RECLINE

The Shower Commode Rise N Recline lifts the user from 55 to 100 cm and has a weight capacity of 325 kg Safe Working Load. It is the world's first shower commode for bariatric users that automatically can reach a 20 degree tilt position by pushing only one button. The chair tilts slowly while lifting the user and reaching a 20 degree tilt in the highest position.

A truly ergonomic multi-purpose hygiene chair

The XXL-Rehab Shower Commode Rise N Recline was developed in response to the growing trend towards showering in bariatric care and skincare as traditionally used methods are uncomfortable for the user and impose strain on the caregiver. The innovative ergonomic design of the shower commode is highly versatile and can also support better routines for washing/showering, hair washing, wound care, toileting etc.



Enhancing lives

The XXL-Rehab Shower Commode Rise N Recline can enhance the quality of life by making assisted showering and other hygiene routines not only safer and more comfortable, but also a more dignified experience.

Making tasks easier and safer

Ergonomic design, battery-powered assistance and the Care Raiser function removes the stress and strain from showering and other everyday hygiene routines, allowing the caregiver to work in a safe and ergonomically sound way.

Improving care, efficiency and communication

This shower commode raises the efficiency of care because it enables a single caregiver to perform the whole showering procedure, as well as a wide range of everyday hygiene tasks. The high level of support provided by the shower commode results in streamlined routines that offer considerable time savings compared with traditional methods. It also enables the caregiver and the user to maintain eye contact for better communication.

- The Care Raiser function is unique, as it lifts the lower part of the user's body to enable easy and discreet access for sensitive hygiene tasks.
- The high backrest and sturdy armrests provide optimum user support and security.
- The reclined position gives full access for the caregiver to most parts of the body and is very comfortable for the user.



Product video
Scan the code and see
the product in use



In lowest position with standard footrests

STRONGLY ENGINEERED AND FUNCTIONALLY DESIGNED EXCLUSIVELY FOR BARIATRIC USERS!

SPECIFICATIONS

Maximum weight capacity	325 kg
Seat depth adjustable	42 / 48 / 54 cm
Total length (lowest position)	140 cm
Total width	74 cm
Seat width	66 cm
Seat height	51-69 cm
Electric tilt	0-20°
Tilt of seat independent from height adjustment	Yes
Width between standard armrests	61 cm
Width between butterfly armrests	76 cm
Width toilet opening	Diameter 24 cm, 35 cm long
Seat height electric adjustable	55-100 cm
Distance from floor to seat underside	44.5 cm
Pushbar from floor	100 cm
Seat to footrest	39-42-45-48 cm
Height seat to armrest	22 cm
Weight	46 kg
Swing-away armrests	Yes
Swing-away footrests	Yes
Removable footrests	Yes
Height adjustable footrests	Yes
Headrest	Yes
4 wheels with individual locks	Yes
Control	Via the control panel
Capacity up/down fully charged battery	Approx. 15 times
Turning radius	156 cm
Frame	Powder lacquered chromium molybdenum steel
Wheels	Diameter 150 mm
Seat / back / armrests	Polyurethane
IP - classification	IP44
ISO - classification	ISO 09 12 03 - 23

ORDERING


Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Shower Commode Rise N Recline incl. standard footrests and charger	66	42 / 48 / 54	0156-066-000
Accessories			
Butterfly armrests (add 14 cm to seat width)			0150-099-012
Headrest adjustable			0156-099-020
Commode seat, 2 part			0150-099-020
Commode bucket			0140-099-001
Footrest right			0150-099-008
Footrest left			0150-099-009
Calf support / amputation support right			0156-099-030
Calf support / amputation support left			0156-099-031



The handset gives the caregiver control over battery-powered functions



Commode bucket and seat



Headrest



Sales • Hire • Service

SHOWER BENCH

STRONGLY ENGINEERED AND FUNCTIONALLY DESIGNED EXCLUSIVELY FOR BARIATRIC USERS!

The armrests are constructed so they go beyond the seat. This makes it easy both to get a good grip when a person wants to sit down but also to get up again. We have designed the armrests to bend out like a butterfly, so the body doesn't get squeezed, getting in and out the shower bench. The butterfly armrests of the XXL-Rehab Shower bench also provide much better functionality for the user, compared to other shower benches with straight and short armrests where it is almost impossible for a bariatric user to get a good grip.

The height of the shower bench is adjustable from 44 to 59 cm. It is very important to adjust the shower bench to the right height to have success getting in and out of the bench. If the shower bench is too high for the user, the feet will lift off the ground, making it hard to get seated in the shower bench. If the shower bench is too low, the user will have trouble getting up.

It is often a problem for bariatric users to be seated comfortably. Sometimes the seat depth is too short, sometimes too deep and with no back support. We have designed the seat to be individually adjusted to provide a comfortable seating. The backrest can be moved +/- 6 cm and the bended tubing assures that the right comfort can be found.

It is easy to clean the XXL-Rehab Shower bench with its removable seat pad. No tools are needed.

It is also made easy to both handle and move the shower bench from place to place. In the construction of the shower bench we have only used high-tech lightweight tubing, to provide the highest strength at the lowest weight.

- Adjustable seat height
- Adjustable seat depth



ORDERING

Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Shower Bench	61	40-52	0141-061-000
XXL-Rehab Shower Bench	71	40-52	0141-071-000

SPECIFICATIONS

Model 1 Model 2

Maximum weight capacity is	325 kg	
Seat width	61 cm	71 cm
Total length	80-101 cm	
Seat depth	40-52 cm	
Seat height	44-59 cm	
Backrest height	37 cm	
Distance from floor to seat underside	33-48 cm	
Total outer width of front legs	69 cm	79 cm
Total height	79-94 cm	
Total weight	15 kg	16 kg
Total outer width on armrests	73 cm	83 cm
Seat / back / armrests	Polyurethan	
Frame	Powder-lacquered chromium molybdenum steel	
ISO - classification	ISO 09 33 03	

BEDSIDE COMMODE

STRONGLY ENGINEERED AND FUNCTIONALLY DESIGNED EXCLUSIVELY FOR BARIATRIC USERS!

The armrests are constructed so they go beyond the seat. This makes it easy both to get a good grip when a person wants to sit down but also to get up again. We have designed the armrests to bend out like a butterfly, so the body doesn't get squeezed, getting in and out the chair. The butterfly armrests of the XXL-Rehab Bedside commode also provide much better functionality for the user, compared to other chairs with straight and short armrests where it is almost impossible for a bariatric user to get a good grip.

The height of the chair is adjustable from 44 to 59 cm. It is very important to adjust the chair to the right height to have success getting in and out of the chair. If the chair is too high for the user, the feet will lift off the ground, making it hard to get seated in the chair. If the chair is too low, the user will have trouble getting up.

It is often a problem for bariatric users to be seated comfortably. Sometimes the seat depth is too short, sometimes too deep and with no back support. We have designed the seat to be individually adjusted to provide a comfortable seating. The backrest can be moved +/- 6 cm and the bended tubing assures that the right comfort can be found.

It is easy to clean the XXL-Rehab Bedside commode with its removable seat pad. No tools are needed and the bedpan is redesigned for a smoothly hygienic cleaning operation.

It is also made easy to both handle and move the chair from place to place. In the construction of the chair we have only used high-tech lightweight tubing, to provide the highest strength at the lowest weight.

- Adjustable seat height
- Adjustable seat depth



ORDERING

Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Bedside commode	61	40-52	0140-061-000
XXL-Rehab Bedside commode	71	40-52	0140-071-000
Accessories			
Commode bucket			0140-099-001

SPECIFICATIONS

Model 1 Model 2

Maximum weight capacity is	325 kg	
Seat width	61 cm	71 cm
Total length	80-101 cm	
Seat depth	40-52 cm	
Seat height	44-59 cm	
Backrest height	37 cm	
Distance from floor to seat underside	33-48 cm	
Distance from floor to bucket underside	21-36 cm	
Total outer width of front legs	69 cm	79 cm
Total height	79-94 cm	
Total weight	14 kg	15 kg
Total outer width on armrests	73 cm	83 cm
Seat / back / armrests	Polyurethan	
Frame	Powder-lacquered chromium molybdenum steel	
ISO - classification	ISO 09 33 03	



MINIMAXX - FOLDING WHEELCHAIR

The XXL-Rehab Minimaxx folds up easily and is easy to bring along. The folding system is very easy to activate, and gives at the same time great support to the seat. The design of the folding system makes the chair rigid – when unfolded the chair is also much more stable than normal foldable chairs.

STW (Stiffness To Weight) was our goal. We used very light tubing technology as used on high end MTB and road bikes. This Minimaxx only weighs approx. 9% of its weight capacity and still provides a smooth drive thanks to the rigid folding system and light tubing.

The XXL-Rehab Minimaxx drives very well due to the position of the rear and front wheels. It is designed so the rear wheels carry most of the user's weight making it much easier for the user to push. The castor wheels are further in front as on standard chairs. This prevents the chair from tipping.

The low seat height makes it easy to access the chair and still propel it.

The push handles are designed so the attendant can easily use his/her body to turn and push the chair. The padded handles are at an angle that puts as little strain to the wrists and shoulders as possible.

The XXL-Rehab Minimaxx has disc brakes as an additional option. Little effort is required to brake the chair even with a full load. To the disc brake is also added a parking brake and the chair also has a wheels lock. Another option is the butterfly armrests. Butterfly armrests add more than 10 cm to the seating width.

The footrests are designed to carry 100 kg each! The footrests are also side/width adjustable. This enables users with legs aligned or with a wide leg spread to use the chair. The very long armrests make it easy for the user to get in and out of the chair.

- Folding rigid frame
- Super easy to push
- Different armrests available



Product video
Scan the code and see
the product in use

STRONGLY ENGINEERED AND FUNCTIONALLY DESIGNED EXCLUSIVELY FOR BARIATRIC USERS!

SPECIFICATIONS

Maximum weight capacity is	325 kg
Seat width	56 / 61 / 66 / 71 cm
Total width	Seat width + 22 cm
Total width folded	35 cm
Total height	103 cm
Total length without footrests	70 cm
Seat depth	46 cm
Seat height	43 cm
Backrest height	51 cm
Weight	Approx. 27 kg
ISO – classification	ISO 12 22 03 – 08

ORDERING

Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Minimaxx Standard incl. footrests	56	46	0100-056-000
XXL-Rehab Minimaxx Standard incl. footrests	61	46	0100-061-000
XXL-Rehab Minimaxx Standard incl. footrests	66	46	0100-066-000
XXL-Rehab Minimaxx Standard incl. footrests	71	46	0100-071-000
XXL-Rehab Minimaxx with disc brakes, incl. footrests	56	46	0109-056-000
XXL-Rehab Minimaxx with disc brakes, incl. footrests	61	46	0109-061-000
XXL-Rehab Minimaxx with disc brakes, incl. footrests	66	46	0109-066-000
XXL-Rehab Minimaxx with disc brakes, incl. footrests	71	46	0109-071-000
Accessories			
Butterfly armrests, pair			0100-099-005
Sport armrests, pair			0100-099-006
Table			0100-099-001
Cane holder for 1 cane			0100-099-022
IV holder and O2 holder, diameter 10 cm			0100-099-020
Calf support / amputation support, right			0100-099-015
Calf support / amputation support, left			0100-099-016
Footrest swing-away, right			0100-099-010
Footrest swing-away, left			0100-099-011




Disc brakes



Table



IV holder and O2 holder



Hand brake operated by attendant



Calf support / amputation support



Sales • Hire • Service

MINIMAXX WITH PUSH MOTOR

The XXL-Rehab Minimaxx With Push Motor is a strong power pack that assists the attendant to push the chair with no effort. It is now easy to move a bariatric patient around. With the XXL-Rehab Minimaxx With Push Motor the caregiver avoids sore shoulders and stress on wrists.

Even on slopping surface up to 5 % the XXL-Rehab Minimaxx With Push Motor proudly helps over you over the top and going down you find the security of our disc braking system. If an immediate stop is needed just pull the lever and the wheelchair stops at once.

The push motor is easy to operate. The control panel is placed on the handle and offers Drive and Reverse mode as well as speed control from 1-5 km/h.

As for the construction of the wheelchair STW (Stiffness To Weight) was our goal. We used very light tubing technology as used on high end MTB and road bikes. This XXL-Rehab Minimaxx only weighs approx. 9% of its weight capacity and still provides a smooth drive thanks to the light tubing.

The castor wheels are further in front as on standard chairs. This prevents the chair from tipping.

The very long armrests make it easy for the user to get in and out of the chair. The low seat height makes it easy to access the chair.

The footrests are designed to carry 100 kg each! The footrests are also side/width adjustable. This enables users with legs aligned or with a wide leg spread to use the chair.

- Push motor is easy to operate
- Different armrests available



STRONGLY ENGINEERED AND FUNCTIONALLY DESIGNED EXCLUSIVELY FOR BARIATRIC USERS!

SPECIFICATIONS

Maximum weight capacity is	325 kg
Seat width	56 / 61 / 66 / 71 cm
Total width	Seat width + 22 cm
Total height	103 cm
Total length without footrests	70 cm
Seat depth	46 cm
Seat height	43 cm
Backrest height	51 cm
Weight	Approx. 27 kg
Weight push motor	23 kg
Weight battery	9 kg
Battery	12V 12AH x 2
Motor	24V 450W 4.5 amps
Tire size	50 mm x 200 mm (Solid Tire) x 2
Charger capacity	Off Board (1.8A)
Range per charge	4.5 km
Max. speed	5 km/h
Gradient	5°
ISO – classification	ISO 12 23 12 - 02

ORDERING

Model	Seat width cm	Seat depth cm	Order No
XXL-Rehab Minimaxx With Push Motor, incl. disc brakes and footrests	56	46	0110-056-000
XXL-Rehab Minimaxx With Push Motor, incl. disc brakes and footrests	61	46	0110-061-000
XXL-Rehab Minimaxx With Push Motor, incl. disc brakes and footrests	66	46	0110-066-000
XXL-Rehab Minimaxx With Push Motor, incl. disc brakes and footrests	71	46	0110-071-000
Accessories			
Butterfly armrests, pair			0100-099-005
Sport armrests, pair			0100-099-006
Table			0100-099-001
Cane holder for 1 cane			0100-099-022
IV holder and O2 holder, diameter 10 cm			0100-099-020
Calf support / amputation support, right			0100-099-015
Calf support / amputation support, left			0100-099-016
Footrest swing-away, right			0100-099-010
Footrest swing-away, left			0100-099-011

Push motor



Disc brakes




Standard armrests



Butterfly armrests



Sport armrests



Sport armrests, side view



Sales • Hire • Service

STAND TALL

The XXL-Rehab Stand Tall is a multi-functional aid; It can be used for walking, mobilising from seated to standing, it is safe to use, gives independence and has an IV option.

The XXL-Rehab Stand Tall has 2 built-in handgrips on the side which allow the user to stand up in small steps from the bed or chair being in full control of the situation. This minimises the risk of falling and therefore also minimises the risk of injury to caregiver or patient himself. Once upright the user can walk around.



The underarm-padding is adjustable in width and angle giving optimal support to the user. This function is supported in full length of the bar.

In addition, it has many handgrip-combinations and easy adjustment, everything to make the XXL-Rehab Stand Tall easy to maneuver.

The brakes are mounted onto the bar, useful in stops of course, but also in speed control and turning: the XXL-Rehab Stand Tall is a safe zone. The brakes can be locked, i.e. when getting up/sitting down.

Using 32 mm tubing ensures an almost forever lasting frame, but also gives a much more simple look and it can be cleaned in few seconds.

IV holder and holder for urine-bag are optional.

The XXL-Rehab Stand Tall is mainly designed for hospital use, but can be used at home as well.

- Height adjustable
- Maximum weight capacity up to 325 kg

SPECIFICATIONS

Maximum weight capacity is	325 kg
Width	76 cm
Adjustable height	100-130 cm
Internal width	71 cm
External width between armrests	55 cm
Depth	72 cm
Front wheel	100 mm
Rear wheel	125 mm
Frame	Powder-lacquered chromium molybdenum steel
ISO - classification	ISO 12 06 12 - 06

ORDERING

Model	Width cm	Length cm	Order No
XXL-Rehab Stand Tall	76	73	0180-076-000
Accessories			
IV holder for Stand Tall			0180-099-008

PATIENT/DINNER CHAIR

STRONGLY ENGINEERED AND FUNCTIONALLY DESIGNED
EXCLUSIVELY FOR BARIATRIC USERS!

The armrests are constructed so they go beyond the seat. This makes it easy both to get a good grip when a person wants to sit down but also to get up again. We have designed the armrests to bend out like a butterfly, so the body doesn't get squeezed, getting in and out the chair. The butterfly armrests of the XXL-Rehab Patient/Dinner chair also provide much better functionality for the user, compared to other chairs with straight and short armrests where it is almost impossible for a bariatric user to get a good grip.

The height of the chair is adjustable from 44 to 59 cm. It is very important to adjust the chair to the right height to have success getting in and out of the chair. If the chair is too high for the user, the feet will lift off the ground, making it hard to get seated in the chair. If the chair is too low, the user will have trouble getting up.

It is often a problem for bariatric users to be seated comfortably. Sometimes the seat depth is too short, sometimes too deep and with no back support. We have designed the seat to be individually adjusted to provide a comfortable seating. The backrest can be moved +/- 6 cm and the bended tubing assures that the right comfort can be found.

It is made easy to both handle and move the chair from place to place. In the construction of the chair we have only used high-tech lightweight tubing, to provide the highest strength at the lowest weight.



SPECIFICATIONS

Model 1 Model 2

Maximum weight capacity is	325 kg	
Seat width	61 cm	71 cm
Total length	80-101 cm	
Seat depth	40-52 cm	
Seat height	44-59 cm	
Backrest height	37 cm	
Distance from floor to seat underside	33-48 cm	
Total width	69 cm	79 cm
Total height	79-94 cm	
Total weight	16.6 kg	17.6 kg
Total outer width on armrests	73 cm	83 cm
Seat / back / armrests	Polyurethan	
Frame	Powder-lacquered chromium molybdenum steel	
ISO - classification	ISO 09 33 03	

ORDERING

Model	Seat width cm	Seat depth cm	Seat height cm	Order No
XXL-Rehab Patient/Dinner Chair	61	40-52	44-59	0160-061-000
XXL-Rehab Patient/Dinner Chair	71	40-52	44-59	0160-071-000



Rehabilitation & Mobility Equipment



Brisbane

10 Devlan St, Mansfield Qld 4122

Phone: 1300 743 710

Gold Coast

7 Dover Dr, Burleigh Heads Qld 4220

Phone: 1300 894 142

Email: **info@qldrehab.com.au**

www.qldrehab.com.au