


NEW



Ref. 1272.HSE06

23G HEAVY SILICONE OIL EXTRACTION CANNULA



Dr. Andreas Mohr,

head of Ophthalmology Department at St. Joseph - Stift Hospital, Bremen, Germany.

“The 10mm cannula is the optimum length to maintain contact with the silicone oil bubble whilst delivering a steady flow of oil during extrusion. Maintaining contact with the bubble is critical to removing heavy oil in a single step and avoid breaks in the bubble leaving smaller oil bubbles on the retina. This cannula is the result of collaboration between DORC and VR surgeons – showing how small product enhancements can make life easier for VR surgeons!”

- **For extraction of heavy silicone oil**
- **10mm long metal tip**
- **With Luer lock connector**