



# PAD A CHEEK<sup>®</sup>

Hose Covers

## 6' Zippered Hose Cover



Instructions on Use and Care  
To install [PAD A CHEEK Hose Cover](#):

Remove the hose from your mask and machine. Unzip the hose cover most of the way. Place your hose in the cover and zip it up. Secure the ends of the hose cover to the hose with the Velcro tabs. Attach your hose to the machine and the mask.

[PAD A CHEEK<sup>®</sup>](#) is made from 100% polyester micro fleece, and can be machine or hand washed and dried (warm water, low heat or air dry). We recommend washing the cover with the zipper closed.

PAD A CHEEK 434-985-4003

157 Sunset Drive, Stanardsville, VA 22973

email [admin@padacheek.com](mailto:admin@padacheek.com) web site [www.padacheek.com](http://www.padacheek.com) Designed and Sewn in the USA

