

PAD A CHEEK[®]

Strap Pad + Dreamwear Top of Head Pad



- Lay the **Top of Head Pad** under the Dreamwear hose connection tube with the pad in the center.
- You can place it with the fleece pad toward your head or away from your head.
- Wrap the Velcro tabs over the Dreamwear and secure.
- Smooth down the **Top of Head Pad** so there are no wrinkles.

PAD A CHEEK Dreamwear Top of Head Pad is 100% polyester microfiber and fleece and can be machine or hand washed.
If you have questions, please contact us.

PAD A CHEEK LLC 434-985-4003 157 Sunset Drive, Stanardsville, VA 22973
Email: admin@padacheek.com web site: www.padacheek.com Designed and Sewn in the USA

