



# PAD A CHEEK<sup>®</sup>

## Strap Pad AirFit P10 with Velcro

- Lay one **Strap Pad** under a side strap with the label touching the inside of the mask strap.
- Wrap the pad around the side and fasten the Velcro together, lining up the edges of the pad so there is no overlap. Make sure that none of the Velcro is sticking out where it can catch your hair or pillowcase or scratch your face when you sleep.
- Repeat for the other strap.

**PAD A CHEEK AirFit P10** is made from 100% polyester micro fleece and can be machine or hand washed in warm water and air or machine dried (low heat).

PAD A CHEEK 434-985-4003  
157 Sunset Drive, Stanardsville, VA 22973

email [admin@padacheek.com](mailto:admin@padacheek.com)

web site [www.padacheek.com](http://www.padacheek.com)

Made in the USA

