



PAD A CHEEK[®]



Stability Strap Mask Stabilizer

- Slide the strap under one of the lower straps on your headgear. Make certain the hook Velcro is on the side away from your face.
- Wrap the strap under your chin and under the other lower strap.
- Secure the Velcro with enough tension to keep the mask from riding up at night.
- You can overlap the strap on the bottom, and attach the hook to the smaller loop section if needed.

PAD A CHEEK Mask Stabilizer is made from 100% polyester micro fiber, and can be machine or hand washed and dried (warm water, low heat or air dry).

PAD A CHEEK 434-985-4003

157 Sunset Drive, Stanardsville, VA 22973

Email: admin@padacheek.com

web site: www.padacheek.com

Designed and sewn in the USA

