



PAD A CHEEK[®]

Ear Loop Protectors Ear Loop Covers



- Wrap the [Ear Loop Covers](#) around the ear loops or nasal Oxygen cannula tubing and secure with the narrow Velcro strips making sure there are not bumps along the edge.
- Make certain that the fold of the fabric is toward your head and the edges of the fabric over your ear.
- Package contains 2 Ear Loop Covers

[PAD A CHEEK Ear Loop Cover](#) is made from 100% polyester micro fleece and can be machine or hand washed with warm soapy water. Do not use bleach or softener.

PAD A CHEEK 434-985-4003

157 Sunset Drive, Stanardsville, VA 22973

Email: admin@padacheek.com website: www.padacheek.com

Designed and manufactured in the USA



On the Swift FX Bella Loops Mask



On the Oxygen Cannula